RIGID TQ CASE FOR SOFTT/SOFTT-W/SAM XT - RIGID SOFTT/SOFTT-W CASE W/MALICE CLIP, BLACK

The RIGID TQ Case for the SOF®TT/SOF®TT-W can hold the SOF®TT or the SOF®TT-Wide tourniquets. The Eleven 10® RIGID TQ Case® is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the SOF®TT/SOF®TT-W utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded Case to stand up to rugged LE/MIL usage Drainage Hole at the bottom of the case Medical Identifying PVC Cross Patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Available Mounting Options - Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt), or MOLLE/PALS Attachment (Tactical Tailor Malice Clip) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10's products are designed, manufactured and assembled in the USA.

Attributes

- Name: RIGID SOFTT/SOFTT-W CASE W/MALICE CLIP, BLACK
- Manufacturer: ELEVEN 10
- Product no.: 100043725
- Mfr. No.: E10-3011M-BLK
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 76mm
- Shipping width: 76mm
- Shipping length: 127mm
- UPC: 700615668243

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RIGID TQ Case Safety Instructions

Introduction

Thank you for choosing the RIGID TQ Case for the SOF®TT/SOF®TTW. This case is designed to securely hold your tourniquet and provide quick access in emergency situations. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Always ensure that the RIGID TQ Case is mounted securely before use.
- Regularly inspect the case for any signs of wear or damage.
- Keep the case and its contents out of reach of children to prevent accidental use.
- Do not modify the case or its components in any way.
- If the case is damaged or does not perform as expected, discontinue use immediately and seek assistance.

Specific Safety Precautions for Use

- Ensure that the tourniquet is properly seated within the case before use.
- Be aware of your surroundings when deploying the tourniquet to avoid injury to yourself or others.
- Practice deploying the tourniquet from the case in a safe environment to familiarize yourself with the process.
- Use the tourniquet only as instructed and for its intended purpose.
- In case of blood exposure, follow appropriate protocols for decontamination and personal safety.

Instructions for Installation and Usage

Installation

1. **Mounting Options:** Choose a suitable mounting option for your RIGID TQ Case. It can be attached to:
 - Police Duty Belt (BladeTech TEKLOK Belt Attachment fits up to a 2.25" belt)
 - MOLLE/PALS Attachment (Tactical Tailor Malice Clip)
 - Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately)
2. **Secure the Hardware:** Use the included Vibratite to secure the mounting hardware in place once you have obtained the desired configuration.
3. **Adjust Height:** Ensure that the case is mounted at a comfortable height for quick access. It can be adjusted to various positions: vertically, horizontally, canted, or inverted.

Usage

1. **Accessing the Tourniquet:** In an emergency, quickly grab the exposed portion of the tourniquet and pull to release it from the case.
2. **Deployment:** Follow the manufacturer's instructions for properly applying the tourniquet to ensure effective use.
3. **PostUse:** After use, clean the case and tourniquet according to the manufacturer's recommendations.

Disposal Instructions

- Dispose of the RIGID TQ Case and its components in accordance with local regulations.
- Do not dispose of the case in regular household waste if it is damaged or poses a safety risk.

Contact Information for Further Support

For any inquiries regarding safety, usage, or support, please refer to the manufacturer's contact details provided with the product packaging.

By following these safety instructions, you can ensure that you use the RIGID TQ Case effectively and safely. Thank you for your attention to these important guidelines.

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