TOURNIQUET/SELF-AID POUCH W/MOLLE ATTACHMENT MULTICAM

The Eleven 10® Tourniquet / Self-Aid Pouch can hold a SOF-T®T Wide Tourniquet and accommodates a QuikClot® Combat Gauze Dressing. It can be worn on your duty belt (vertical or horizontal) or you can utilize the included Malice Clips to mount it (vertical only) to any MOLLE/PALS webbing. The pouch has a cross patch to identify it as "medical" and allows you to carry the life-saving equipment you need on your belt where it belongs. Holds the SOF-T®T Wide Tourniquet with or without QuikClot® Combat Gauze Medical Identifying PVC Cross Patch included Velcro® and Snap Closure Fits a Duty Belt or can be attached to MOLLE/PALS using the included Malice Clips Can be mounted vertically or horizontally on belt, vertical only with Malice Clips Made with 1000D Cordura All of Eleven 10®'s products are designed, manufactured and assembled in the USA.

Attributes

- Name: TOURNIQUET/SELF-AID POUCH W/MOLLE ATTACHMENT MULTICAM
- Manufacturer: ELEVEN 10
- Product no.: 100043742
- Mfr. No.: E10-1005M-MTC
- Size: 1 Person
- Delivery weight: 0.136kg
- Shipping height: 51mm
- Shipping width: 127mm
- Shipping length: 178mm
- UPC: 700615669769

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Safety Instruction Guide for TOURNIQUET/SELF AID POUCH

Introduction

Thank you for choosing the Eleven 10® Tourniquet / SelfAid Pouch with MOLLE Attachment. This guide provides essential safety information and instructions to ensure the safe use and maintenance of your product. Please read this guide carefully before using the pouch.

General Safety Guidelines

- Ensure the pouch is kept in a dry and clean environment to maintain the integrity of the materials.
- Regularly inspect the pouch for any signs of wear or damage. Do not use if damaged.
- Keep out of reach of children unless under adult supervision.
- Always store the tourniquet and other medical supplies in a manner that prevents unauthorized access.
- In case of emergency, ensure that all users are familiar with how to access and use the pouch.

Specific Safety Precautions for Use

- The pouch is designed to hold a SOFT®T Wide Tourniquet and QuikClot® Combat Gauze Dressing. Ensure that only compatible items are stored inside.
- When wearing the pouch on your duty belt, ensure that it is securely fastened to prevent accidental removal or loss.
- If using Malice Clips for MOLLE/PALS attachment, ensure they are properly secured and that the pouch is mounted in a manner that does not obstruct movement.
- Do not overstuff the pouch, as this may compromise its structural integrity and effectiveness.
- In case of a medical emergency, follow established medical protocols for the use of tourniquets and other medical supplies.

Instructions for Installation and Usage

Installation

1. Duty Belt Attachment:

- Position the pouch on your duty belt.
- It can be worn either vertically or horizontally based on your preference.
- Ensure the pouch is securely fastened and does not slide or shift.

2. MOLLE/PALS Attachment:

- Use the included Malice Clips for vertical mounting only.
- Align the clips with the MOLLE/PALS webbing and thread them through.
- Pull firmly to secure the pouch in place.

Usage

1. Accessing the Tourniquet:

- Open the pouch using the Velcro® and snap closure.
- Remove the SOFT®T Wide Tourniquet and any additional contents as needed.

2. Applying the Tourniquet:

- Follow standard first aid procedures for applying the tourniquet.
- Ensure that the tourniquet is applied above the injury site and tightened appropriately.

Disposal Instructions

- Dispose of the pouch and its contents in accordance with local regulations for medical waste.
- If the pouch is damaged, do not attempt to repair it. Instead, dispose of it responsibly.
- Always ensure that any medical supplies are disposed of in a safe manner to prevent accidental exposure.

Contact Information for Further Support

For any inquiries regarding the safety or use of the Eleven 10® Tourniquet / SelfAid Pouch, please refer to the contact information provided with your product.

Thank you for your attention to these safety guidelines. Your safety and the effective use of this product are our top priorities.

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