PTAKS RUGGED REFILL POD - PTAKS RUGGED-REFILL POD W/COMPRESSED GAUZE BLACK

Contents (1) CLEER Medical Trauma Dressing 4" (1) Phokus Research Frog Gauze OR QuikClot Combat Gauze LE (1) Nitrile Gloves (1) Phokus Research Group / Eleven 10 Frog Tape (1) North American Rescue Hyfin Vent Compact Chest Seal, Twin Pack Pouch Sold Seperate

Attributes

- Name: PTAKS RUGGED-REFILL POD W/COMPRESSED GAUZE BLACK
- Manufacturer: ELEVEN 10
- Product no.: 100043751
- Mfr. No.: PTAKSGABLK
- Size: 1 Person
- Delivery weight: 0.245kg
- Shipping height: 197mm
- Shipping width: 38mm
- Shipping length: 178mm
- UPC: 700615672325

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PTAKS RUGGED REFILL POD Safety Instruction Guide

Introduction

Thank you for choosing the PTAKS RUGGED REFILL POD. This product is designed to provide essential first aid supplies in emergency situations. It is important to follow the safety guidelines outlined in this manual to ensure safe and effective use.

General Safety Guidelines

- Always keep the PTAKS RUGGED REFILL POD out of reach of children.
- Inspect the contents regularly to ensure they are intact and not expired.
- Use the product only for its intended purpose, as described in the instructions.
- In case of an emergency, seek professional medical help immediately after using the first aid supplies.
- Report any unsafe products or incidents to the appropriate authorities.

Specific Safety Precautions for Use

- Ensure that all items are stored in a cool, dry place to maintain their effectiveness.
- Avoid using any supplies if the packaging is damaged or if the items appear to be contaminated.
- Wear nitrile gloves when handling any of the supplies to reduce the risk of infection.
- Follow the instructions provided with each component for proper application and usage.
- Be aware of any allergies or sensitivities to materials included in the kit.

Instructions for Installation and Usage

1. Preparation:

- Open the pouch carefully to avoid damaging the contents.
- Lay out all items on a clean surface.

2. Using the Trauma Dressing:

- Remove the dressing from its packaging.
- Place the dressing over the wound and apply gentle pressure.
- Secure the dressing in place as needed.

3. Using the Gauze:

- Choose either the Phokus Research Frog Gauze or QuikClot Combat Gauze LE.
- Pack the gauze into the wound, ensuring it is applied firmly.
- If using QuikClot, follow the specific instructions on the packaging for application.

4. Applying the Chest Seal:

- Remove the chest seal from its packaging.
- Apply it over any open chest wound, ensuring a tight seal to prevent air from entering the chest cavity.

5. Using Frog Tape:

- Use the Frog Tape to secure dressings or gauze in place as necessary.
- Ensure the tape does not cover any breathing passages.

6. Gloves:

- Always wear nitrile gloves when handling any medical supplies to prevent contamination.

Disposal Instructions

- Dispose of used or expired medical supplies in accordance with local regulations.
- Do not dispose of medical waste in regular trash; use designated medical waste disposal services if available.
- Ensure that any sharp objects (e.g., used needles) are disposed of in a proper sharps container.

Contact Information for Further Support

For any inquiries regarding the PTAKS RUGGED REFILL POD or to report safety concerns, please reach out to your local distributor or regulatory authority.

By following these guidelines, you can ensure the safe and effective use of the PTAKS RUGGED REFILL POD. Thank you for prioritizing safety and preparedness.

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