SABA POUCH RANGER GREEN

The Eleven 10® SABA (Self-Aid / Buddy-Aid) Pouch is considered the “upgraded version” of our Tourniquet / Self-Aid Pouch and can hold a SOF-T®T Wide Tourniquet, C-A-T® Tourniquet, or a SWAT-TTM. It can be worn on your duty belt or you can utilize the included Malice Clips to mount it to any MOLLE/PALS webbing. The pouch has a cross patch to identify it as “medical” and allows you to carry the life-saving equipment you need on your belt where it belongs. There are 2 snap positions which adjust the height of the pouch to accommodate different tourniquet sizes. There is a small pocket inside the flap which can hold a pair of exam gloves. Features: Holds the SOF-T®T Wide Tourniquet, C-A-T® Tourniquet (Gen 7 and Older Generations), or a SWAT-TTM Medical Identifying PVC Cross Patch Included Velcro® and Snap Closure Fits a Duty Belt or can be Attached to MOLLE/PALS Using the Included Malice Clips Made with 500D Cordura All of Eleven 10®'s products are designed, manufactured and assembled in the USA.

Attributes

- Name: SABA POUCH RANGER GREEN
- Manufacturer: ELEVEN 10
- Product no.: 100043762
- Mfr. No.: E10-1019-RGR
- Size: 1 Person
- Delivery weight: 0.136kg
- Shipping height: 51mm
- Shipping width: 127mm
- Shipping length: 178mm
- UPC: 700615672097

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Safety Instructions for the Eleven 10® SABA Pouch

Introduction

Thank you for choosing the Eleven 10® SABA Pouch. This pouch is designed to safely carry essential medical equipment, including various types of tourniquets. To ensure safe and effective use, please read and follow these safety instructions carefully.

General Safety Guidelines

- Ensure the pouch is used only for its intended purpose – to carry medical equipment.
- Inspect the pouch regularly for any signs of wear or damage. Do not use if damaged.
- Keep the pouch out of reach of children and vulnerable individuals unless supervised.
- Store the pouch in a cool, dry place to maintain the integrity of the materials.
- Follow all local regulations and guidelines regarding the use of medical equipment.

Specific Safety Precautions for Use

- Always ensure that the pouch is securely closed when in use to prevent accidental opening.
- Familiarize yourself with the contents of the pouch and their proper usage before an emergency arises.
- Avoid overloading the pouch beyond its intended capacity to ensure ease of access to medical supplies.
- Use caution when handling sharp or potentially hazardous items stored within the pouch.
- If using the pouch in a highstress environment, ensure it is easily accessible and identifiable.

Instructions for Installation and Usage

1. Mounting the Pouch:

- The SABA Pouch can be attached to a duty belt or mounted on MOLLE/PALS webbing using the included Malice Clips.
- To attach to a duty belt:
 - Slide the pouch onto the belt until it fits snugly.
- To mount on MOLLE/PALS:
 - Thread the Malice Clips through the webbing and secure them firmly.

2. Adjusting the Pouch Height:

- The pouch features two snap positions to adjust its height.
- To adjust:
 - Unsnap the pouch from its current position.
 - Move the pouch to the desired height and snap it securely in place.

3. Using the Pouch:

- Open the flap to access the contents.
- Use the small pocket inside the flap to store a pair of exam gloves or similar items.
- Ensure that all items are securely placed in the pouch before closing.

Disposal Instructions

- Dispose of the pouch in accordance with local regulations regarding textile and medical waste.
- If the pouch is damaged or no longer needed, ensure it is disposed of in a way that prevents accidental access to its contents.

Contact Information for Further Support

For any inquiries or concerns regarding the Eleven 10® SABA Pouch, please refer to the appropriate contact points provided by your retailer or manufacturer.

By following these safety instructions, you can help ensure the safe and effective use of the Eleven 10® SABA Pouch. Thank you for prioritizing safety in your medical preparedness.

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