ULTRA EAR BTE RECHARGEABLE- SINGLE - ULTRA EAR BTE RECHARGEABLE HEARING AID 2PK

Rechargeable 48 Hours battery life on a single charge Lightweight, compact behind-the-ear design can be used in either ear On/Off Switch dB Limiting @ 105dB Can be used in either ear Includes 3 sizes foam ear tip 3 bottom cell batteries Brush and wax pick / Screwdriver USB charging cable. Rechargeable 48 Hours battery life on a single charge Lightweight, compact behind-the-ear design can be used in either ear On/Off Switch dB Limiting @ 105dB Noise Reduction Rating 22dB Can be used in either ear 3 sizes foam ear tip 3 bottom cell batteries Brush and wax pick / Screwdriver USB charging cable.

Attributes

- Name: ULTRA EAR BTE RECHARGEABLE HEARING AID 2PK
- Manufacturer: WALKERS GAME EAR
- Product no.: 430100937
- Mfr. No.: GWP-RCHUE-2PK
- Color: Black
- NRR: -
- Style: -
- Delivery weight: 0.077kg
- Shipping height: 34mm
- Shipping width: 127mm
- Shipping length: 203mm
- UPC: 888151023372

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ULTRA EAR BTE RECHARGEABLE SAFETY INSTRUCTIONS

Introduction

Thank you for choosing the ULTRA EAR BTE RECHARGEABLE hearing aid. This guide provides essential safety instructions and information to ensure safe and effective use of your hearing aid. Please read this guide carefully before using the product.

General Safety Guidelines

To ensure your safety and the proper functioning of your hearing aid, please adhere to the following guidelines:

- Always use the product as directed.
- Keep the hearing aid out of reach of children and pets.
- Avoid exposing the hearing aid to moisture, extreme temperatures, or direct sunlight.
- Do not attempt to disassemble or modify the hearing aid.
- Regularly check for any signs of damage or wear and replace parts as necessary.
- Dispose of batteries and the hearing aid in accordance with local regulations.
- Report any unsafe conditions or accidents to the appropriate authorities.

Specific Safety Precautions for Use

When using the ULTRA EAR BTE RECHARGEABLE hearing aid, please observe the following precautions:

- Ensure the hearing aid is charged before use to avoid sudden power loss.
- Use only the provided USB charging cable for recharging.
- Do not use the hearing aid if it is damaged or not functioning correctly.
- Be cautious when adjusting the volume to avoid potential hearing damage.
- If you experience discomfort or pain while using the hearing aid, discontinue use immediately and consult a professional.
- Keep foam ear tips clean and replace them as needed to maintain hygiene.

Instructions for Installation and Usage

Follow these instructions to ensure proper installation and usage of your hearing aid:

1. Charging the Hearing Aid:

- Connect the USB charging cable to a power source.
- Plug the other end into the charging port of the hearing aid.
- Allow the hearing aid to charge fully (up to 48 hours of battery life on a single charge).

2. Inserting the Hearing Aid:

- Choose the appropriate foam ear tip size for comfort.
- Gently insert the hearing aid behind your ear, ensuring a snug fit.
- Adjust the ear tip to ensure it is comfortable and secure in your ear.

3. Turning On/Off:

- To turn on, slide the On/Off switch to the "On" position.
- To turn off, slide the switch to the "Off" position.

4. Adjusting the Volume:

- Gradually increase the volume to a comfortable level using the volume control.
- Avoid sudden increases in volume to protect your hearing.

5. Maintenance:

- Clean the hearing aid regularly with the provided brush and wax pick.
- Store the hearing aid in a dry place when not in use.

Disposal Instructions

Proper disposal of the hearing aid and its components is essential for environmental safety:

- Remove and dispose of batteries according to local regulations.
- Do not throw the hearing aid in regular household waste.
- Check for any recycling programs in your area for electronic devices.

Contact Information for Further Support

For any inquiries regarding safety or product performance, please reach out to the designated contact point provided by the manufacturer or retailer.

Thank you for your attention to these safety instructions. Enjoy your experience with the ULTRA EAR BTE RECHARGEABLE hearing aid!

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